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
大井 悦雅 e\_ohi00@shimadzu-techno.co.jp

## Chlordecone (Kepone)

This analytical method produces for the determination of chlordecone in breast milk, maternal blood and cord blood by liquid-chromatography/tandem quadrupole mass spectrometry (LC/MS/MS).

[Breast milk]

Five grams of breast milk sample spiked 2 ng of  $^{13}\text{C}_{10}$ -chlordecone as surrogate is added 10 mL of diethyl ether and 2 mL of saturated potassium oxalate aqueous solution. And then, extract with 20 mL of ethanol/n-hexane (1:1) by liquid-liquid extraction for 30 minutes. This process is repeated with 20 mL of n-hexane. All extracts are collected and washed with 20 mL of 5% sodium chloride solution and dehydrate with sodium sulfate anhydrous. A rotary evaporator concentrated the extract until dryness. After measurement of fat content, residue is dissolved with n-hexane to 10 mL. Two milliliter of raw extract is treated with acetonitrile/n-hexane participation. The purified extract is concentrated to dryness under a gentle nitrogen stream and dissolved the residue with 1 mL of methanol. If analytes have opacity, filtrate with PTFE filter. The analytes are determined in the multiple-reaction-monitoring (MRM) mode as the precursor/product ion pair of  $m/z$  506.7/426.5 for chlordecone and  $m/z$  516.8/435.7 for  $^{13}\text{C}_{12}$ -chlordecone. The method detection limit (MDL) and the method quantification limit (MQL) are 3.1 and 7.9 pg/g, respectively. The average of recoveries ( $N=6$ ) from 0.5 ng chlordecone added breast milk is 98.8%, and the relative standard division is 4.0%.



[Maternal blood]

Twenty grams of maternal blood sample spiked 2 ng of  $^{13}\text{C}_{10}$ -chlordecone as surrogate is added 12 mL of saturated ammonium sulfate aqueous solution. And then, extract with 48 mL of ethanol/n-hexane (1:3) by liquid-liquid extraction for 30 minutes. This process is repeated twice with 40 mL of n-hexane. All extracts are collected and washed with 40 mL of n-hexane-washed water and dehydrated with sodium sulfate anhydrous. A rotary evaporator concentrated the extract solution until dryness. After measurement of fat content, residue is dissolved with n-hexane to 10 mL. Two milliliter of raw extract is treated with acetonitrile/n-hexane participation. The purified extract is concentrated to dryness under a gentle nitrogen stream and dissolved the residue with 1 mL of methanol. If analytes have opacity, filtrate with PTFE filter. The analytes are determined in the multiple-reaction-monitoring (MRM) mode as the precursor/product ion pair of  $m/z$  506.7/426.5 for chlordecone and  $m/z$  516.8/435.7 for  $^{13}\text{C}_{12}$ -chlordecone. The method detection limit (MDL) and the method quantification limit (MQL) are 0.81 and 2.1 pg/g, respectively. The average of recoveries (N=6) from 0.5 ng chlordecone added control serum is 103%, and the relative standard division is 7.1%.

[Cord blood]

Thirty grams of cord blood sample spiked 2 ng of  $^{13}\text{C}_{10}$ -chlordecone as surrogate is added 18 mL of saturated ammonium sulfate aqueous solution. And then, extract with 72 mL of ethanol/n-hexane (1:3) by liquid-liquid extraction for 30 minutes. This process is repeated twice with 60 mL of n-hexane. All extracts are collected and washed with 60 mL of n-hexane-washed water and dehydrated with sodium sulfate anhydrous. A rotary evaporator concentrated the extract solution until dryness. After measurement of fat content, residue is dissolved with n-hexane to 10 mL. Two milliliter of raw extract is treated with acetonitrile/n-hexane participation. The purified extract is concentrated to dryness under a gentle nitrogen stream and dissolved the residue with 1 mL of methanol. If analytes have opacity, filtrate with PTFE filter. The analytes are determined in the multiple-reaction-monitoring (MRM) mode as the precursor/product ion pair of  $m/z$  506.7/426.5 for chlordecone and  $m/z$  516.8/435.7 for  $^{13}\text{C}_{12}$ -chlordecone. The method detection limit (MDL) and the method quantification limit (MQL) are 0.54 and 1.4 pg/g, respectively.

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