2013.8.12

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The screenshot shows the 'ComSpec exe for Japanese OS (64 bit version)' window. It contains several sections for file selection and configuration. Five numbered steps are overlaid on the left side of the window:

- ① Select .cdf file for ComSpec processing. (File... button)
- ② Select database .csv file for ComSpec processing. (File... button)
- ③ Select save folder and input file name. Chromat pictures (.jpg) and a result file (.csv) are created. (File... button)
- ④ - R Package instillation -  
☐ Install Packages. Check that you are online!  
☒ Package ncodf are installed
- ⑤ - Check your R version -  
R-3.0.0 Input your R version for using.  
☒ Place of R folder is [Program files]. Check off if C drive directly

Below these steps, there are tabs for 'Separation parameter' and 'Mass parameter'. The 'Separation parameter' tab is active, showing a dropdown menu set to 'Two-dimensional' with the text 'Select dimension.' below it. Further down, there are input fields for 'MPeriod: The modulation time period.' (value: 4) and 'Phaseshift: The time of phase shift.' (value: -1.3). At the bottom right, there is a 'ComSpec run' button.

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The screenshot shows a software window with two tabs: "Separation parameter" and "Mass parameter". The "Separation parameter" tab is active. Inside this tab, there is a dropdown menu currently set to "Two-dimensional" with a downward arrow, followed by the text "Select dimension.". Below this, there are two input fields. The first is labeled "MPeriod: The modulation time period." and contains the value "4". The second is labeled "Phase.shift: The time of phase shift." and contains the value "-1.3". At the bottom of the window, there is a button labeled "ComSpec run".

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The screenshot shows the 'Mass parameter' tab of the ComSpec software. It contains several input fields and a dropdown menu for configuring mass spectrometry parameters. The 'Accurate MS' dropdown is selected. The 'CountMS' field is set to 5. The 'Tolerance (m/z)' is calculated as  $\text{HR.wide.order} \times \text{HR.wide.order.fluc}$ . The 'HR.wide.order' field is set to 0.01, and the 'HR.wide.order.fluc' field is set to 5. The 'MSthreshold' field is set to 0. A 'ComSpec run' button is located at the bottom right of the interface.

Separation parameter   Mass parameter

Accurate MS   Select MS type.

5   CountMS: The number of MS counting for searching/assignment of target compounds.

Tolerance (m/z) = HR.wide.order × HR.wide.order.fluc

0.01   HR.wide.order: The decimal point position in accurate mass.

5   HR.wide.order.fluc: The degree of fluctuation.

0   MSthreshold: The minimum threshold of ion intensity.

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