

Instructions of a 'MDF LM'

Shunji Hashimoto

National Institute for Environmental Studies,

Tsukuba, JAPAN, 305-8506

Tel&fax: +81-29-850-2531

E-mail: shunji@nies.go.jp

Read me first

This program was developed for study purposes, so please be forewarned that you are using it at your own risk. We are not able to provide compensation for any damage it might cause.

The program requires a Microsoft Windows 10 or later “64 bit” environment. It also needs installing a netCDF4 (64bit) run time liberally. The liberally can be obtained from unidata web site (<https://www.unidata.ucar.edu/software/netcdf/>).

BE SURE TO BACK UP YOUR INPUT DATA BEFORE USING THIS PROGRAM.

Please let us know if you have any special requirements or suggestions. We will consider them for inclusion in our development process.



This work is licensed under a [Creative Commons Attribution-NonCommercial 4.0 International License](https://creativecommons.org/licenses/by-nc/4.0/).

Our contact information;

National Institute for Environmental Studies, Onogawa 16-2, Tsukuba, JAPAN, 305-8506


Tel&fax: +81-29-850-2531

e-mail: shunji@nies.go.jp

Reference: Hashimoto *et al.*(2013) Journal of Chromatography A, 1282, 183-189

What is a “MDF LM”?

This software can...



Interface of a “MDF LM”

- ◆ directly read a netCDF file, which is a common format for MS data,
- ◆ and extract mass spectra in range of mass defects,
- ◆ under various mass accuracy (mass resolution) conditions,
- ◆ and finally, output the result as netCDF.

Main settings for data extraction

1) **Input File Name***: choose filename(s) from the dialog box that will appear after right-clicking.


3) **Threshold**: check and input an integer from 0 to 1,000,000 if needed. Recommended values are 500, 1,000 or 3,000 depending on the situation. (Commas are not required.)

2) **Mass Defect Filter***: The value must be from 0 to 1. (See **Technical Guide 1.**)


The screenshot shows the 'MDF LM v1.0.0' dialog box. It has a 'Close' button at the top left. Below it is the 'Input File Name:' label followed by a large text input field and a browse button (three dots). Below this are two rows of settings: the first row has a checkbox for 'Threshold' (unchecked) and a text input field with '0'; the second row has a text input field for 'Mass Defect Filter' with '0' and a dropdown for 'Mass Accuracy (u)' set to '0.03'. A note '(absolute defect value, from 0 to 1)' is next to the mass accuracy dropdown. At the bottom right are 'Execute' and 'Cancel' buttons. At the very bottom, there are labels for 'Input File Name:', 'Lines:', and 'Time(min):'.

4) **Mass Accuracy***: choose a value from 1 to 0.001u. (This parameter will determine the accuracy of the mass, e.g., MA: 0.05u is equivalent to a mass resolution of $m/\Delta m$ 6,000 @ m/z 300.)


*all fields followed by an asterisk must be filled in.



After all settings are complete, click the **Execute** button.




Running program
Elapsed time is displayed.





Completed process

If the process has successfully completed, an output file, formatted as netCDF, is created in the same folder as the input file.

Output file: The file name is automatically created by appending the parameters used for the processing.



 080508-09-CENTROIDED.cdf_STCorrected 3-10Cl 0.1u TH_500 MD_on IE_10 MR_FullRange VIP_on NLS_off.cdf	908,473 KB
 080508-09-CENTROIDED.cdf_STCorrected.cdf	908,473 KB

Input file: The file must be formatted in netCDF.

1: Filter using mass deficiency

Mass excess (Mass defect):

The mass excess of a nuclide is the difference between its **actual mass** and its **mass number** in atomic mass units.

— from Wikipedia

Mass deficiency:

Mass deficiency is a state where the mass excess is negative. (**actual mass** < **mass number**) N and atoms whose mass number exceeds ¹²C exhibit mass deficiency .

Mass deficiency is used for simple data filtration before Cl or Br isotopic pattern checking.

Specifically, mass spectra showing mass deficiency within a 0 to -0.2 range are left in the data and the other spectra are removed from the data.

	¹ H	⁴ He	...	¹² C	¹⁴ N	¹⁶ O	...	³⁵ Cl	...	⁷⁹ Br
Atom No.	1	4		12	14	16		35		79
Exact mass	1.007825	4.0026		12	14.003074	15.9949		34.96885268		78.9183371
defference	+0.007825	+0026		±0	+0.003074	-0.0051		-0.03114732		-0.0816629

February 8, 2019