
2020.06.12

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The screenshot shows the NMFdeconvolution exe (64 bit version) window. It has a title bar with standard Windows window controls. The interface is divided into several sections:

- Choose files -**: Contains a label "Select .cdf file for NMFdeconvolution" and a "File..." button. This is annotated with a circled 1.
- Select save files position -**: Contains a label "Select save folder and input file name. Chromat pictures (.jpg) and a result file (.csv) are created." and a "File..." button. This is annotated with a circled 2.
- Required R Packages -**: Contains text: "Check that all of packages [ncdf4, EBImage, xtable, NMF] are installed. If not installed, run the enclosed script in the software folder." This is annotated with a circled 3.
- Check your R version -**: Contains a text input field with "R-4.0.0" and a label "Input your R 4version for using." This is annotated with a circled 4.
- NMF parameter**: A tabbed interface with "Data properties and Prefilters" selected. It contains:
  - A dropdown menu for "Select NMF algorithm." with "Frobenius" selected. This is annotated with a circled 5.
  - A dropdown menu for "Select initial seeding method." with "nndsvd" selected.
  - A text input field for "Factor setting: The number of factors (ranks)." with "5" entered.
  - A text input field for "Output setting: The number of factors (ranks) to output. Should be lower than the Factor setting" with "4" entered.
  - A text input field for "Precision: Precision of m/z value in NMFdeconvolution" with "0.1" entered.
  - A text input field for "Peak picking parameter: Watershed method is applied. The highest resolution is 1 (Integer)" with "1" entered.
- Buttons**: At the bottom, there is an "NMFdeconvolution run" button and a checkbox labeled "Execute database construction". The "run" button is annotated with a circled 6.

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The screenshot shows the 'NMF parameter' tab of the *I\_NMFdeconvolution\_expand.exe* application. The interface is divided into two main sections: 'NMF parameter' and 'Data properties and Prefilters'. The 'NMF parameter' section contains the following settings:

- ⑤-1** Select NMF algorithm: Frobenius (dropdown menu)
- ⑤-2** Select initial seeding method: nndsvd (dropdown menu)
- ⑤-3** Factor setting: The number of factors (ranks): 5 (input field)
- ⑤-4** Output setting: The number of factors (ranks) to output. Should be lower than the Factor setting: 4 (input field)
- Precision: Precision of m/z value in NMFdeconvolution: 0.1 (input field)
- Peak picking parameter: Watershed method is applied. The highest resolution is 1 (Integer): 1 (input field)

At the bottom of the interface, there is a button labeled 'NMFdeconvolution run' and a checkbox labeled 'Execute database construction'.

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The screenshot shows the 'NMF parameter' tab of the software interface. It contains several input fields and checkboxes for configuring the deconvolution process. Numbered callouts are placed over the interface to guide the user:


- ⑤-5** points to the 'MS type' dropdown menu, which is currently set to 'HRscan'.
- ⑤-6** points to the 'MDF value' input field, which is set to '0.2'.
- ⑤-7** points to the 'MStreshold' input field, which is set to '0'.
- ⑤-8** points to the 'Range of m/z' input fields, which are set to '30' and '600'.
- ⑤-9** points to the 'MPeriod' input field, which is set to '4'.
- ⑥** points to the 'Execute database construction' checkbox, which is checked.

At the bottom of the interface, there is a button labeled 'NMFdeconvolution run' and a checkbox labeled 'Execute database construction'.

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
GCxGC Mixture Touch:

[http://www.mixture-platform.net/Mixture\\_Touch\\_open/](http://www.mixture-platform.net/Mixture_Touch_open/)

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The screenshot shows the 'Library Search submission' window. On the left, a vertical line with four numbered boxes (0, 1, 2, 3) indicates the sequence of steps. Step 0 points to a caution message. Step 1 points to the 'File...' button. Step 2 points to the 'Expected calculation time (h)' input field. Step 3 points to the 'Library Search run' button.

Library Search submission

0 Caution!! This program overwrite the files of NIST MS Search (AUTOIMP.MSD, secondLocator.fil) to default. Save and store them outside of the NIST folder, if you have your own setting.  
- Choose files -

1 Select .txt file for Library Search  
File...

2 Expected calculation time (h) 4.5  
\*Output file will be copied in the folder of input file, after all the search done. Thus, expected time should be longer than NIST processing time. 14000 list (=Max) takes around 4.5 hours. FinalCheck\_forNISTsearch file is preferable to be directly imported in NIST MSsearch, manually.

3 C:/NIST17 Input folder path of NISTlibrary including correct version of NIST to call up  
Library Search run

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```
Unknown: Scan 1 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 616; RMF: 648; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 2 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 619; RMF: 652; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 3 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 617; RMF: 650; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 4 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 623; RMF: 656; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 5 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 613; RMF: 646; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 6 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 623; RMF: 657; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 7 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 620; RMF: 653; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 8 0      Compound in Library Factor = N/A
```

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The screenshot shows the DBCreator+ application window (64 bit version) with the following sections and steps:

- Step 1:** Select .txt file of NIST MS search output. A "File..." button is highlighted.
- Step 2:** Select .csv file of MSpeaklist result. A "File..." button is highlighted.
- Step 3:** Select save folder and input file name. Chromat pictures (.jpg) and a result file (.csv) are created. A "File..." button is highlighted.
- Step 4:** Check your R version - R-4.0.0 Input your R version for using. ☒ Place of R folder is [Program files]. Check off if C drive directly
- Step 5:** NIST MS Search setting GC setting
  - Match Factor (MF) threshold: NIST Hit list with the setting MF value and over is extracted. Range 0 ~ 1000. Value: 800
  - Number of hits to print: Should be same with NIST MS search setting. Value: 1
  - Keyword: Keyword in chemical formula to extract from the hit list. e.g.) Chlorinated compound is extracted by Cl

At the bottom right, there is a button labeled "DBCreator+ run".

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NIST MS Search setting GC setting

⑤-1 900 Match Factor (MF) threshold: NIST Hit list with the setting MF value and over is extracted. Range 0 ~ 1000

⑤-2 1 Number of hits to print: Should be same with NIST MS search setting

⑤-3 Keyword: Keyword in chemical formula to extract from the hit list. e.g.) Chlorinated compound is extracted by Cl

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NIST MS Search setting GC setting

⑤-4 1 GC1 tolerance (min): Multiply entries in the hit list within this range is removed.

1 GC2 tolerance (sec): Same as above. Both GC1 and GC2 tolerance is considered at the same time.

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