


12th June, 2020



Operation manual for GUI-based NMFwithDBcreator

Developer

National Institute for Environmental Studies

Guest Scientist

Yasuyuki ZUSHI

Framework for NMFwithDBcreator

- NMFwithDBcreator involves a four step process and consists of the following three executable files:
 - *I_NMFdeconvolution_expand.exe*
 - *II_IV_LibrarySearch.exe*
 - *III_DBCreator+.exe*
- Each of the files is executed in the corresponding step; then, the database file is created as the final output.
- Detailed explanations are given on the following pages.

Requirements for NMFwithDBcreator

- OS: Windows 10
- At least 10 GB RAM and a 64-bit OS are recommended for large amounts of data.
- Installation of “R” statistical software is required. Freely downloadable from:
<https://cran.ism.ac.jp/>
- netCDF data from GCxGC-HRTOFMS have been confirmed to work in this software.
- Do not touch the computer keyboard during a few seconds, because this tool automatically inserts the required information – there is no need for you to do anything. **Please save and close important files before starting the tool.**
- **Demonstration data are included in this software folder. Also, it is available via the following link (Japanese site with English translation):**
<https://www.nies.go.jp/analysis/downloads.html>
- NIST MS Search software is required for “*II_IV_LibrarySearch.exe*” and “*III_DBCreator+.exe*”.

I_NMFdeconvolution_expand.exe

Step 1:

Spectral deconvolution is executed. Then, all the generated peaks, including the original peak with RTs and MS spectra, are listed in a csv file (optional). In addition, a file with the format for a library batch search in the NIST library is generated.

I_NMFdeconvolution_expand.exe

How to Execute (1)

The screenshot shows the NMFdeconvolution exe (64 bit version) window. It has a light gray background and a title bar with standard Windows window controls. The window is divided into several sections. On the left side, there are six numbered circles (1 through 6) pointing to specific parts of the interface. 1. Points to the 'Choose files' section, specifically the 'Select .cdf file for NMFdeconvolution' label and the 'File...' button. 2. Points to the 'Select save files position' section, specifically the 'Select save folder and input file name. Chromat pictures (.jpg) and a result file (.csv) are created.' label and the 'File...' button. 3. Points to the 'Required R Packages' section, specifically the text 'Check that all of packages [ncdf4, EBImage, xtable, NMF] are installed. If not installed, run the enclosed script in the software folder.' 4. Points to the 'Check your R version' section, specifically the 'R-4.0.0' input field and the 'Place of R folder is [Program files]. Check off if C drive directly' checkbox. 5. Points to the 'NMF parameter' section, specifically the 'Frobenius' dropdown menu for 'Select NMF algorithm.' and the 'hndsvd' dropdown menu for 'Select initial seeding method.' 6. Points to the 'NMFdeconvolution run' button at the bottom right of the window. The 'Data properties and Prefilters' tab is also visible, showing various settings like 'Factor setting: The number of factors (ranks)' set to 5, 'Output setting: The number of factors (ranks) to output' set to 4, 'Precision: Precision of m/z value in NMFdeconvolution' set to 0.1, and 'Peak picking parameter: Watershed method is applied. The highest resolution is 1 (Integer)' set to 1. There is also an unchecked checkbox for 'Execute database construction'.

NMFdeconvolution exe (64 bit version)

- Choose files -

① Select .cdf file for NMFdeconvolution
File...

- Select save files position -

② Select save folder and input file name. Chromat pictures (.jpg) and a result file (.csv) are created.
File...

- Required R Packages -
Check that all of packages [ncdf4, EBImage, xtable, NMF] are installed. If not installed, run the enclosed script in the software folder.

③

- Check your R version -
R-4.0.0 Input your R 4version for using.
☒ Place of R folder is [Program files]. Check off if C drive directly

④

NMF parameter Data properties and Prefilters

⑤ Frobenius Select NMF algorithm.
hndsvd Select initial seeding method.
5 Factor setting: The number of factors (ranks).
4 Output setting: The number of factors (ranks) to output. Should be lower than the Factor setting
0.1 Precision: Precision of m/z value in NMFdeconvolution
1 Peak picking parameter: Watershed method is applied. The highest resolution is 1 (Integer)

⑥

NMFdeconvolution run ☐ Execute database construction

① Choose your measurement data file (.cdf).

② Set an output folder and file name.

③ You need to install the required program packages. If not, install them by running a script in the provided file “Run_me_for_package_installation.r”.

④ Set file path to run “R”. You need to know where and which version of “R” is installed on your computer. Is it located directly on the C drive or in “Program Files”?

⑤ and ⑥ See next two pages.

I_NMFdeconvolution_expand.exe

How to Execute (2)

The screenshot shows the 'NMF parameter' tab of the 'I_NMFdeconvolution_expand.exe' application. The interface includes two tabs: 'NMF parameter' and 'Data properties and Prefilters'. The 'NMF parameter' tab is active, displaying several settings:

- ⑤-1** Select NMF algorithm. The dropdown menu is set to 'Frobenius'.
- ⑤-2** Select initial seeding method. The dropdown menu is set to 'nndsvd'.
- ⑤-3** Factor setting: The number of factors (ranks). The input field is set to '5'.
- ⑤-4** Output setting: The number of factors (ranks) to output. Should be lower than the Factor setting. The input field is set to '4'.
- Precision: Precision of m/z value in NMFdeconvolution. The input field is set to '0.1'.
- Peak picking parameter: Watershed method is applied. The highest resolution is 1 (Integer). The input field is set to '1'.

At the bottom of the window, there is a button labeled 'NMFdeconvolution run' and a checkbox labeled 'Execute database construction' which is currently unchecked.

- ⑤-1 Select NMF algorithm and initial seeding method. Frobenius and nndsvd method are the default settings.
- ⑤-2 Input the number of factors (ranks). “Output setting” describes the number of factors to output. Factors with a higher rank are chosen in descending order.
- ⑤-3 Set precision of m/z, used as the variance in NMF deconvolution: m/z = 100.1, 100.2, 100.3, ... in the case of 0.1. Lower values have higher calculation cost.
- ⑤-4 Peak range is defined based on this parameter. Basically the value one should be chosen. If you want to expand the range for peaks, choose a larger value (must be an integer).

I_NMFdeconvolution_expand.exe


How to Execute (3)

The screenshot shows the 'Data properties and Prefilters' tab of the *I_NMFdeconvolution_expand.exe* application. The interface is divided into two main sections: 'NMF parameter' and 'Data properties and Prefilters'. The 'Data properties and Prefilters' section contains the following elements:

- ⑤-5** Select MS type: A dropdown menu with 'HRscan' selected.
- ⑤-6** Mass Defect filter (MDF): Only for HRscan mode. Radio buttons for 'MDF on' and 'MDF off' (selected).
- ⑤-7** MDF value: The value for MDF (Doesn't work in case MDF off is chosen). Input field with '0.2'.
- ⑤-8** MSthreshold: The minimum threshold of ion intensity. Input field with '0'.
- ⑤-9** Range of m/z: Should be within the range in measurement data. Input fields with '30' and '600'.
- ⑤-9** MPeriod: The modulation time period. Input field with '4'.
- ⑥** A button labeled 'NMFdeconvolution run'.
- ⑥** A checkbox labeled 'Execute database construction' which is checked.

- ⑤-5 Select “HRscan” for data with high mass resolution. Select “scan” for nominal mass data. This parameter is related to MDF only. Usually, “HRscan” should be selected.
- ⑤-6 MDF setting as a pre-filter for the data. The next column describes the MDF value. See [Hashimoto et al., 2013, J. Chromatogr. A 183-189](#) for details of MDF processing.
- ⑤-7 Threshold for ion intensity values is adjusted here.
- ⑤-8 Describes the range of m/z in NMF deconvolution. This range should be within the range of the measurement data.
- ⑤-9 Describes the modulation period (sec) in the measurement by GCxGC.
- ⑥ If checked, files of peak lists and NIST Search are generated in database construction.

Screenshot of Process



- After clicking button “NMFdeconvolution run” with correct configuration, “R” is activated and starts to load the code according to the configuration.
- After loading all the code, peak deconvolution for each peak begins.
- Typically, it takes 1 to 5 seconds for a peak, and an hour for an entire chromatogram of size 100 MB to be deconvoluted.
- The calculation cost increases according to the precision of m/z . It takes an hour to generate three deconvoluted layers with precision $m/z = 1$, an hour to generate a layer with precision $m/z = 0.1$, and 2~3 hours to generate a layer with precision $m/z = 0.05$.
- A CDF file is generated in the save folder as the output of the deconvolution.
- In the following process, files of peak lists and NIST Search are generated if ⑥ is checked.

Output File (1)


The generated CDF file can be opened with the software GCImage.

Other several free tools, such as R, is available.

Easily, the following site provides a web-application to visualize GCxGC data.

GCxGC Mixture Touch:

http://www.mixture-platform.net/Mixture_Touch_open/



NMF-based deconvolution

$$Y \approx WH$$
$$\min_{W, H \geq 0} [D(Y, WH) + R(W, H)]$$

The layer is output according to the “output setting”.

Output File (2)

- A number of output files are generated depending on the setting for “output setting”.
- The name of the original cdf file is set to *“filename_layer0.cdf”*.
- The name of the deconvoluted file for the first layer is set to *“filename_layer1.cdf”*.
- If ⑥ is checked, files of peak lists and NIST Search are generated following the peak picking process, after the deconvolution process. This process takes a few minutes.
- Both *“MSpeaklist_filename_layer0.csv”* and *“MSpeaklist_forNISTsearch_filename_layer0.txt”* are generated in the process.
- Furthermore, *“Combined_MSpeaklist_filename.csv”* and *“Combined_MSpeaklist_forNISTsearch_filename.txt”* form the combined list of information on peaks in each layer.

II_IV_LibrarySearch.exe

Step 2 and Step 4:

The file for batch search in NIST library is loaded and the result is output.

Step 4 is provided merely as a final check of extracted peaks after the database construction in step 3.

* Step 4 using this exe file is not absolutely essential. The simplest way is to manually open and load the file for batch search in the NIST MS Search software.

II_IV_LibrarySearch.exe

How to Execute (1)

The screenshot shows a window titled "Library Search submission" with the following content:

- 0** Caution!! This program overwrite the files of NIST MS Search (AUTOIMP.MSD, secondLocator.fil) to default. Save and store them outside of the NIST folder, if you have your own setting.
- Choose files -
- 1** Select .txt file for Library Search
File...
Expected calculation time (h) 4.5
*Output file will be copied in the folder of input file, after all the seach done. Thus, expected time should be longer than NIST processing time. 14000 list (=Max) takes around 4.5 hours. FinalCheck_forNISTsearch file is preferable to be directly imported in NIST MSsearch, manually.
- 2** C:/NIST17 Input folder path of NISTlibrary including correct version of NIST to call up
- 3** Library Search run

- 0 The following two files are automatically placed in the folder for the NIST library program, prior to starting the search:
"AUTOIMP.MSD" and "secondLocator.fil"
- **Caution!!** Original files are stored in the location "C:/NIST_{version}/MSSEARCH"; please move the original files to a different location, if you have your own search setting.
- ① Select a file to execute batch search in the NIST library. The file name should be "Combined_MSpeaklist_forNISTsearch_filename.txt", if the user has not changed it.

II_IV_LibrarySearch.exe

How to Execute (2)

The screenshot shows a window titled "Library Search submission" with the following content:

- 0** Caution!! This program overwrite the files of NIST MS Search (AUTOIMP.MSD, secondLocator.fil) to default. Save and store them outside of the NIST folder, if you have your own setting.
- Choose files -
- 1** Select .txt file for Library Search
File...
Expected calculation time (h) 4.5
*Output file will be copied in the folder of input file, after all the seach done. Thus, expected time should be longer than NIST processing time. 14000 list (=Max) takes around 4.5 hours. FinalCheck_forNISTsearch file is preferable to be directly imported in NIST MSsearch, manually.
- 2** C:/NIST17 Input folder path of NISTlibrary including correct version of NIST to call up
- 3** Library Search run

- ② Input an estimated required time to complete NIST Search. Owing to the process of automatically moving an output file after completion of NIST Search, the required time should be estimated with plenty of margin; e.g., it generally takes 4 hours for 14000 entries, and thus 4.5 hours is recommended as input.
- ③ Input the folder path of the NIST library to use in the search. The path depends on which version of the NIST library is used.
- All parameter settings in the NIST library should be completed before execution. For the parameter settings, run the NIST MS Search software and select "library search option". The recommended settings are "similarity" with "simple search" in the search mode, and one for "Number of hits to print".

Output File

```
Unknown: Scan 1 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 616; RMF: 648; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 2 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 619; RMF: 652; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 3 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 617; RMF: 650; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 4 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 623; RMF: 656; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 5 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 613; RMF: 646; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 6 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 623; RMF: 657; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 7 0      Compound in Library Factor = N/A
Hit 1 : <<Hexasiloxane, 1,1,3,3,5,5,7,7,9,9,11,11-dodecamethyl->>;<<C12H38O5Si6>>; MF: 620; RMF: 653; Prob: -1.00; CAS:995-82-4; Mw: 430; Lib: <<mainlib>>; Id: 39868.
Unknown: Scan 8 0      Compound in Library Factor = N/A
```

An example of the output

- A text file describing the results of the NIST Search of all peaks is the output.
- Search results with a number matching the setting “Number of hits to print” are assigned to a peak.
- In the case of the upper figure, “Number of hits to print” is set as one, and therefore, the search result is recorded in the output file.
- The output file name is set to:
“Combined_MSpeaklist_forNISTsearch_filename.txtOutputLibrarySearch.txt”.

III_DBCreator+.exe

Step 3:

The list of all peaks obtained by deconvolution is extracted under certain conditions in this step.

Through this process, the database file and the file for the final check of spectra in the NIST library are generated simultaneously.

III_DBCreator+.exe

How to Execute (1)

The screenshot shows the DBCreator+ application window with the following sections and steps:

- Step 1:** Select .txt file of NIST MS search output. A "File..." button is highlighted.
- Step 2:** Select .csv file of MSpeaklist result. A "File..." button is highlighted.
- Step 3:** Select save folder and input file name. Chromat pictures (.jpg) and a result file (.csv) are created. A "File..." button is highlighted.
- Step 4:** Check your R version - R-4.0.0 Input your R version for using. ☒ Place of R folder is [Program files]. Check off if C drive directly
- Step 5:** NIST MS Search setting GC setting
 - Match Factor (MF) threshold: NIST Hit list with the setting MF value and over is extracted. Range 0 ~ 1000. Value: 800
 - Number of hits to print: Should be same with NIST MS search setting. Value: 1
 - Keyword: Keyword in chemical formula to extract from the hit list. e.g.) Chlorinated compound is extracted by Cl

At the bottom right, there is a button labeled "DBCreator+ run".

① Select a txt file for the NIST search result. The name of the file is "Combined_MSpeaklist_forNISTs earch_filename.txtOutputLibraryS earch.txt", if no changes have been made to the file name in the previous step.

② Select a csv file with the peak list; default name is "Combined_MSpeaklist_filename.csv".

③ Choose save file location and input save file name.

④ Select version and location of software R.

III_DBCreator+.exe

How to Execute (2)

NIST MS Search setting GC setting

⑤-1 900 Match Factor (MF) threshold: NIST Hit list with the setting MF value and over is extracted. Range 0 ~ 1000

⑤-2 1 Number of hits to print: Should be same with NIST MS search setting

⑤-3 Keyword: Keyword in chemical formula to extract from the hit list. e.g.) Chlorinated compound is extracted by Cl

DBcreator+ run

- ⑤-1 Lower threshold value of MF to extract from the peak list.
- ⑤-2 “Number of hits to print”, which is selected in “Library Search Option” in NIST MS Search.
- ⑤-3 Entries containing this keyword in the assigned formula are extracted from the list. If “Cl” is input, entries with “Cl” in their formula and at the same time, greater than the MF threshold, are extracted.

III_DBCreator+.exe

How to Execute (3)

The screenshot shows the 'GC setting' tab of the 'NIST MS Search setting' window. It contains three input fields with their respective labels and descriptions:

- ⑤-4** Input field: 1. Label: GC1 tolerance (min): Multiply entries in the hit list within this range is removed.
- ⑤-4** Input field: 1. Label: GC2 tolerance (sec): Same as above. Both GC1 and GC2 tolerance is considered at the same time.
- ⑤-5** Input field: 4. Label: MPeriod: The modulation time period.

At the bottom of the window is a button labeled 'DBcreator+ run'.

- ⑤-4 Input tolerance of retention time in GC1 and GC2.
- If 1 min and 1 sec are chosen, entries in the same assignment within this tolerance are regarded as duplicates. Thus, the entries are deleted except for the entry with the highest MF.
- ⑤-5 Input modulation time period applied in the measurement.

Output File

BlobID	Compound	Group	Rtc	Inte	Peak I (mir	Peak II (se	IS	分子式	MF	RMF	CAS	LibID	Lib	file.n	deconv.nur	MS1	MS2
1	Heptadecan	NA	0	0	10.88133	1.230071	0	C18H38	917	934	13287-23-	22596	mainlib	fileN	0	57.07812	71.09436
NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	10.8538	10.77631
2	Hexadecan	NA	0	0	18.948	0.753914	0	C20H40O2	962	966	111-06-8	20871	mainlib	fileN	0	56.06208	57.06959
NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	12.67844	7.087579
3	Octadecan	NA	0	0	23.348	0.119039	0	C22H44O2	916	916	123-95-5	20912	mainlib	fileN	0	56.06212	57.06972
NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA	12.02491	7.305549

Part of the list in database file

Information on a peak is summarized across two rows

...

The following three output files are generated.

1. The list of extracted peaks based on certain conditions. This file consists of meta data such as compound name, RT, formula. The file name is set to "SimplePeaklist_conditionKeyword_MF_MFvaluefilename.csv".
2. This file comprises a compound database that is used by the TSEN program. It includes MS spectra in addition to file (1), and the file name is "Database_conditionKeyword_MF_MFvaluefilename.csv".
3. This file displays the compound spectra in the final list of the database. The file name is set to "FinalCheck_forNISTsearchDatabase_conditionKeyword_MF_MFvaluefilename.csv.txt". The MS spectra in the final list are visually confirmed point-by-point by loading this file in NIST MS Search.

Precautions

- About the source file
 - Keep the source files “1NMFdeconvolution.r”, “2MSpeaklist.r”, “3FinalDBcreate.r” and “4DBfinalCheckInNIST.r” in the same directory as the .exe files “I_NMFdeconvolution_expand_en.exe”, “II_IV_LibrarySearch.exe”, and “III_DBCreator+_en.exe”.
- If the software freezes, you can exit by clicking the icon at the bottom right of the screen.
- License: Artistic License 2.0
- Disclaimer: The developer is not responsible for any hardware or software damages, data losses or false inferences caused by the use of this software.