2013.8.12

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The screenshot shows the T-SEN exe special edition (64 bit version) window. It has a blue title bar and a light beige background. The interface is divided into several sections with numbered steps on the left:

- ①** - Choose files -
Select .cdf file for T-SEN processing
- ②** Select database .csv file for T-SEN processing
- ③** - Select save files position -
Select save folder and input file name. Two chromat pictures (.jpg) and a result file (.csv) are created.
- ④** - T-SEN Package installation -
☐ Install Packages. Check that you are online!
☒ Package T-SEN and ncdf are installed
- ⑤** - Check your R version -
R-2.14.2 Input your R version for using.
☐ Place of R folder is [Program files]. Check off if C drive directly
☐ R is 32 bit version. T-SEN with spetial edition is needed.

Below these steps are two tabs: "Separation parameter" (selected) and "Mass parameter". Under "Separation parameter", there is a dropdown menu set to "Two-dimensional" with the text "Select dimension." below it. Below the dropdown are several input fields with labels:

- BaseIS.search: Search range of base IS for RT correction in entire chromatogram (±sec).
- Range.correct: The degree of peak range correction after peak finding (±sec).
- search.range.x-axis: Search range in first GC (±sec). Rounded to MPereod each in 2D.
- Search.range.y-axis: Search range of a peak in second GC (±sec). 2D only
- MPeriod: The modulation time period.
- Phase.shift: The time of phase shift.

At the bottom right of the window is a button labeled "T-SEN run".

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☐ R is 32 bit version. T-SEN with special edition is needed.

Separation parameter Mass parameter

 Select dimension.

 BaseIS.search: Search range of base IS for RT correction in entire chromatogram (\pm sec).

 Range.correct: The degree of peak range correction after peak finding (\pm sec).

 search.range.x-axis: Search range in first GC (\pm sec). Rounded to MPeriod each in 2D.

 Search.range.y-axis: Search range of a peak in second GC (\pm sec). 2D only

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✓ R is 32 bit version. T-SEN with special edition is needed.

Separation parameter **Mass parameter**

Accurate MS Select MS type.

5 CountMS: The number of MS counting for searching/assignment of target compounds.

Tolerance (m/z) = HRwide.order × HRwide.order.fluc

0.01 HRwide.order: The decimal point position in accurate mass.

5 HRwide.order.fluc: The degree of fluctuation.

0 MSthreshold: The minimum threshold of ion intensity.

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